



GC7021-CAUDALIMETRO MANOMETRO

MOD.G108

Caudalímetro con manómetro para baja presión. Suministro de gas regulable entre 0-15 l/min. Salida con pitón para mascarilla o para humidificador

Para uso medicinal.

CARACTERISTICAS TECNICAS

- Clase: IIA
- Normativa:
UNE-EN ISO 15002
- Gas: Gases medicinales
- Conexión de entrada:
AFNOR/SEO/DIN/FRO/CM/BS
- Conexión de salida: 1/4M -
9/16M
- Presión de entrada: 4.5 bar
- Caudal estandar: 0-15 l/min
- Dimensiones: 95 x 90 x 40 mm
- Peso: 325g
- Condiciones ambientales de
utilización: 0°C +40°C
- Condiciones ambientales de
almacenamiento:
-20°C +60°C

ACCESORIOS

- Humidificadores
- Derivadores/Duplicadores
Accesorio selector
- Manguera y pinza de conexión
a barra técnica

CAUDALIMETROS FLOWMETERS



GC7021-FLOWMETER WITH PRESSURE GAUGE

MOD. G108

Flowmeter with pressure gauge for low pressure. Adjustable gas supply between 0-15 l/min. Outlet with hose connector or connector for humidifier.

For medical use.

TECHNICAL CHARACTERISTICS

- Class: IIA
- Standard:
UNE-EN ISO 15002
- Gas: Medical gases
- Inlet connection:
AFNOR/SEO/DIN/FRO/CM/BS
- Outlet connection: 1/4M -
9/16M
- Inlet pressure: 4.5 bar
- Standard flow: 0-15 l/min
- Dimensions: 95 x 90 x 40 mm
- Weight: 325g
- Operating temperature:
0°C +40°C
- Storage temperature:
-20°C +60°C

ACCESSORIES

- Humidifier
- Duplicator accessory
- Selector accessory
- Hose and clamp for connection
to technical bar